



## manugraphy®

*Accurate & reliable hand analysis*

**manugraphy®** enables the analysis of the hand function at highest quality level.

Scan the **pressure distribution** to get a reliable and accurate **analysis of the hand.**


# Hand pressure platform

### manugraphy® key features:

- collect pressure and force data during static & dynamic grasping and holding
- work with reliable, individually calibrated, capacitive sensors to evaluate anatomical segments of the hand individually
- follow pre-defined protocols and create a report within seconds, automatically
- quickly integrate the sensor in your lab or medical environment and sync with other systems



## Technical information



## manugraphy® software features



novel GmbH (Global, GER)  
Ismaninger Str. 51, 81675 Munich  
tel: +49 (89) 417767-0  
e-mail: sales@novel.de  
web: www.novel.de  
copyright © novel GmbH - May 2022

novel electronics inc. (North America)  
964 Grand Avenue St. Paul, MN 55105  
tel: +1 (651) 221-0505  
e-mail: novelinc@novelusa.com  
web: www.novelusa.com